- (6) Plot (D_n² D_m²) with (n-m) and get the slope.
- (7) Also Plot D_n² with n and compare the result.

OBSERVATIONS:

(a) Spherometer readings (Calculation for R)

Pitch of the Spherometer = ... cm

Least Count of the Spherometer = cm

Table 1: Determination of R of the lens

		Spherometer reading		-	7	1 19 10 9 -	
SN		convex surface (a)	Plan surface (b)	$h = a \sim b$	1	$R=f^2/6h + h/2$	
1			N. 4 1.55		1. 22	The state of	
2.	- 0	ALLEY OF	100 m	5 4 . 0	I LA		
3.		The Property of		2	100	1 1 1	79
4.	-	CORN TOWN	16	10300 B	134		14
5.	10	THE PLANT		Tay 19	12.36		THE STATE OF

l is the distance between any two outer legs of the spherometer

(b) Microscope reading

Vernier Constant of microscope = 0.01. cm

R=120cm = 1200 mm

Table 2: Determination of the diameter of rings.

	Microscope reading (Left)			Microscope reading (Right)			Diameter	
SN of rings	M.S.	V.S.	S. MS.+ (V.S.xV.C.)		V.S.	MS.+ (V.S.xV.C.)	$D = a \sim b$	
4 th	47	67	47.65	43	83	43.83	3.86	
19	48	45	48.45	43	00	43	5.45	
12 th	48	91	48.91	42	43	42.43	6.48	
16	49	67	49.67	42	02	42.02	7.65	
20 th	49	99	49.99	41	61	41.61	8.38	

DATA ANALYSIS:

Slope of the plot = $\frac{3.75}{1.00}$, λ (from graph) = $\frac{.7812.5}{1.00}$

Error calculation

SN	n-m	Dn ² · Dm ²	λ	2	$\delta \lambda_i = \lambda_1 - \overline{\lambda}$	δλ, 2	σχ
1.	20-16	11.7019	0.0006094		1.1598x 10	1-34x10-8	2 8
2.	20 - 12	28 - 234	0.00073526	chall Jr 30	-7.18×10-6	5.15×1011	
3	20 8	40. 194	0.00070350	The same	2-1575 XIS	4-78×10+0	
4	20 4	55.632	0.00072437	Carl Tale	1.005 × 10-6	1.01×10-12	
5	16-12	16.5321	0.00086104	7.2538	-1.3566×15-4	1-840x 10-8	1.2355
6	16 - 8	28. 4921	000075052	X 10-4	-2.514×10-5		
7	16-4	43-9301		1 4 4	-3.722×165	1	
8	12 - 8	11.96	00006229	1	1-0.247 15	1.05 ×10-8	
9	12-4	27.398	00007134	8	1.19 × 10-5	1.4161x1010	
10	8-4	15.436	0.00078		-5.462×165	2.983×10 3	
		Σλ = 7, 2	53896 X10	Par dia u	$\Sigma\delta\lambda_i^2 = 1.3$	74×10-6	

RESULTS & CONCLUSION

- The wave length of sodium light = $\bar{\lambda} \pm \sigma_{\lambda} = 7.2538 \times 10^{-4} \pm 2.2355 \times 10^{-4}$
- (v) Standard value of $\lambda = 58.93 \text{ Å} = 5.893 \times 10^4 \text{ m/m}$
- % error = 23.00%.

DISCUSSION

In lab we used the apparatus to find Newton's ring we used the travelling microscope to see the rings and used main and arcular scale for determining the diameter

- Instruments must be handled with care-
- (II) Experiment should be done in donk room.
- (III) The light must fall sufficiently on the microscope

ANSWER THE FOLLOWING QUESTIONS:

A. Before performing Experiment.

(1) What are Newton's rings? How are they formed?

=> Newton's sings are the phenomena in which an interference pattern is created by the reflection of light between two

they are formed by the interference between light waves reflected from top bottom surface of air film form between less and glass-sheet.

(2) Give working formula to light.

= The working formula to light is;

$$\lambda = \frac{Dn^2 - Dm^2}{4(n-m)R}$$

(3) What is interference? What are essential conditions to produce interference

fringes? = The net effect of the combination of two or more waves moving on intersecting or coincident paths.
Conditions to produce interference fringes are:

(a) Source should be coherent and mono-chromatic.

(b) Source should be of same intensity.

(c) Distance between screen should be small.

(4) How circular interference fringes are formed?

= This is due to use of plano-convex lens in Newton's ring.
experiment, Newton's rings are circular because the air film
in between the glass and lens at the point of contact are of equal thickness about a point in a circular way.

(5) What will happen to the rings if the glass plate of the combination is replaced by a plane mirror?

= If the glass plate of the combination is replaced by a plane mirror then, the interference fringes will not be visible as plane missor reflects light uniformly.

B: After Performing Experiment

(1) What would be your observation in transmitted light?

ing appears to be bright.

(2) Why is the central spot dark in your experiment?

=> It is because the light suffers a phase of inverse upon reflecting from a medium with high refraction index.

(3) Why should you use a convex lens of large radius of curvature? Why does the fringes get closer and thinner as we move away from the centre?

= we should use convex lens of large radius of curvature so that when it is placed in our film, the observes will have then a large diameter and the error in the measurement will be less.

The fringes get closer and thinner as we move away from center as moving from center increases the air gap which lowers the radius thus making them closer and thinner.

Date:

Name of teacher:

Teacher's Signature

